

## Co-Dx REGULATORY BULLETIN

**Salt Lake City, September 19<sup>th</sup>, 2022**

**Re:** Response to request for information on the impact of new SARS-CoV-2 variants of interest and variants of concern on the performance of SARS-CoV-2 molecular IVD test kits designed, developed, and manufactured by Co-Diagnostics, Inc. (Salt Lake City, USA)

The World Health Organization and other agencies around the world are currently tracking the predominant circulating and emerging SARS-CoV-2 variants to proactively anticipate their possible impact on public health.

**Table 1** shows a current list of the Variants of Concern (VOC) or Interest (VOI) identified by WHO, MHRA, ECDC, or CDC.

**Table 1 SARS-CoV-2 Genomic Variants Currently Classified as Variants of Concern or Interest**

Variant of Concern / Interest Status by Agency	WHO Nomenclature	Pango Lineage	Nextstrain Clade	Designation, Date	First detected, Date	BEI Reference Material Number
VOI [MHRA]	Delta	B.1.617.2	21A, 21I, 21J	WHO: VOC 11MAY-21 MHRA: VOC APR-2021 ECDC: VOC 24MAY-2021	India, Oct-2020	NR-55611
VOC [WHO, MHRA, ECDC, CDC]	Omicron	B.1.1.529, B.1.1.529.XE, BA.1, BA.1.1, BA.2, BA.3, BA.4, BA.5	21K, 21L, 21M	WHO: VOC 26NOV-2021 MHRA: VOC NOV-2021 ECDC: VOC 26NOV-2021 CDC: VOC 30NOV-2021	Multiple countries, Nov-2021	NR-56496 NR-56520 NR-56803 NR-58616
VOI [MHRA, ECDC]	Omicron	BA.2 + L452X, BA.2.75	22D	MHRA: VOI JUL-2022 ECDC: VOI SEP-2022	Multiple countries, May-2022, JUN-2022	N/A

*Note: VOC – Variant of Concern, VOI – Variant of Interest, WHO – World Health Organization, MHRA – Medicines and health Products Regulatory Agency (UK), ECDC – European Centre for Disease Prevention and Control, CDC – Centers for Disease Control and Prevention*

The WHO has reported that PCR testing continues to effectively detect the B.1.1.529 (Omicron) variant, except for S gene target failure (SGTF) in some PCR assays. (WHO, 2022) SGTF in Omicron is not a concern for Co-Diagnostics assays, which instead target the RdRp and E genes. The WHO also reports that studies are currently ongoing to assess the impact of Omicron on antigen-based rapid diagnostic tests (Ag-RDTs).

Since the initial design of CoPrimers targeting Gene RdRp and Gene E in Feb-2020, Co-Diagnostics has conducted monthly BLASTn queries of subsampled sequences in the Nextstrain database to monitor their homology against emerging strains. These in silico analyses confirm that none of the mutations related to the lineages B.1.1.529, B.1.1.529.XE, BA.1, BA1.1, BA.2, BA.3, BA.4 and BA.5 (Omicron), or Omicron sub-variants BA.2 + L452X or BA.2.75 occur in regions of RdRp or Gene E targeted by any COVID-19 tests designed and manufactured by Co-Diagnostics.

In April 2021, the WHO identified lineage B.1.167.2 (Delta) (PANGO lineages, 2021), which remains a Variant of Interest by designation of MHRA. This Delta lineage contains one point mutation in the region targeted by the RdRp Fwd (Forward) CoPrimer. As shown in **Table 2**, this CoPrimer is used in all COVID-19 tests designed and manufactured by Co-Diagnostics.

**Table 2 COVID-19 Tests Designed and Manufactured by Co-Diagnostics, Inc.**

Test Name	Product Code	Gene RdRp (CoPrimer RdRp Fwd / CoPrimer RdRp Rev)	Gene E (CoPrimer E Fwd / CoPrimer E Rev)
Logix Smart™ Coronavirus Disease 2019 (COVID-19)	COVID-K-001	Yes	No
Logix Smart SARS-CoV-2 (genes RdRp/E)	COVID-K-002	Yes	Yes
Logix Smart ABC (Influenza A/B, SARS-CoV-2)	ABC-K-001	Yes	Yes
	ABC-K-002	Yes	Yes
Logix Smart SARS-CoV-2 DS	COVDS-K-003	Yes	Yes
	COVDS-K-004	Yes	Yes

Since June 2021, Delta expanded to represent a growing proportion of subsampled sequences in the Nextstrain database, peaked in December 2021, and has since diminished to less than 10% of circulating sequences as Omicron has become the dominant circulating strain. This rise and fall in Delta is reflected in the monthly homology statistics for RdRp and E Gene Targets of Co-Diagnostics COVID-19 tests shown in **Table 3** and **Table 4**.

**Table 3 In Silico Analysis Performed Over Time for RdRp Gene Target**

<b>Date of CoDx's Analysis for RdRp Marker</b>	<b>SARS-CoV-2 (sequences in analyzed subsample)</b>	<b>Sequences in the pool with 100% homology</b>	<b>Single nucleotide mutation events: Sequences with 1 mismatch on CoDx target (98% homology)</b>	<b>Double nucleotide mutation events: Sequences with 2+ mismatches on CoDx target (95% homology)</b>	<b>Multiple nucleotide mutation events: Sequences with 3+ mismatches on CoDx target (&lt;95% homology)</b>
27-Jan-20	14	14 (100%)	0 (0%)	0 (0%)	0 (0%)
04-Feb-20	53	53 (100%)	0 (0%)	0 (0%)	0 (0%)
17-Mar-20	571	570 (99.8%)	1 (0.2%)	0 (0%)	0 (0%)
06-Apr-20	3639	3634 (99.86%)	5 (0.14%)	0 (0%)	0 (0%)
04-May-20	4468	4459 (99.80%)	9 (0.2%)	0 (0%)	0 (0%)
03-Jun-20	4558	4537 (99.54%)	21 (0.46%)	0 (0%)	0 (0%)
06-Jul-20	11361	11328 (99.71%)	33 (0.29%)	0 (0%)	0 (0%)
10-Aug-20	22054	22012 (99.81%)	42 (0.19%)	0 (0%)	0 (0%)
09-Sep-20	4417	4394 (99.48%)	23 (0.52%)	0 (0%)	0 (0%)
12-Oct-20	5139	5114 (99.51%)	25 (0.49%)	0 (0%)	0 (0%)
05-Nov-20	3494	3463 (99.11%)	31 (0.89%)	0 (0%)	0 (0%)
04-Dec-20	3407	3371 (98.94%)	36 (1.06%)	0 (0%)	0 (0%)
04-Jan-21	3540	3505 (99.01%)	35 (0.99%)	0 (0%)	0 (0%)
04-Feb-21	3962	3909 (98.66%)	53 (1.34%)	0 (0%)	0 (0%)
01-Mar-21	4024	3960 (98.41%)	64 (1.59%)	0 (0%)	0 (0%)
07-Apr-21	4025	3962 (98.43%)	63 (1.57%)	0 (0%)	0 (0%)
06-May-21	3923	3847 (98.06%)	76 (1.94%)	0 (0%)	0 (0%)
01-Jun-21	3883	3000 (97.86%)	83 (2.14%)	0 (0%)	0 (0%)
06-Jul-21	3883	3566 (91.84%)	317 (8.17%)	0 (0%)	0 (0%)
02-Aug-21	3782	3172 (83.87%)	610 (16.13%)	0 (0%)	0 (0%)
01-Sep-21	3534	2894 (81.89%)	637 (18.02%)	3 (0.08%)	0 (0%)
01-Oct-21	3559	2882 (80.98%)	672 (18.88%)	5 (0.14%)	0 (0%)
01-Nov-21	3572	1648 (46.14%)	1919 (53.72%)	5 (0.14%)	0 (0%)
06-Dec-21	3386	1239 (36.59%)	2140 (63.20%)	7 (0.21%)	0 (0%)
03-Jan-22	3472	1293 (37.24%)	2171 (62.53%)	8 (0.23%)	0 (0%)
01-Feb-22	3201	1361 (42.52%)	1836 (57.36%)	4 (0.12%)	0 (0%)
01-Mar-22	3247	1831 (56.39%)	1411 (43.46%)	5 (0.15%)	0 (0%)
01-Apr-22	2967	1837 (62.01%)	1127 (37.89%)	3 (0.10%)	0 (0%)
01-May-22	2822	1938 (69.04%)	874 (30.61%)	10 (0.35%)	0 (0%)
01-Jun-22	2905	2202 (75.8%)	699 (24.06%)	4 (0.14%)	0 (0%)
01-Jul-22	2817	2570 (91.23%)	245 (8.70%)	2 (0.07%)	0 (0%)
01-Aug-22	2817	2570 (91.23%)	245 (8.70%)	2 (0.07%)	0 (0%)
01-Sep-22	2818	2611 (92.66%)	205 (7.27%)	2 (0.07%)	0 (0%)

**Table 4 In Silico Analysis Performed Over Time for Gene E Target**

Date of CoDx's Analysis for Gene E Marker	SARS-CoV-2 (sequences in analyzed subsample)	Sequences in the pool with 100% homology	Single nucleotide mutation events: Sequences with 1 mismatch on CoDx target (98% homology)	Double nucleotide mutation events: Sequences with 2+ mismatches on CoDx target (95% homology)	Multiple nucleotide mutation events: Sequences with 3+ mismatches on CoDx target <95% homology)
27-Jan-20	14	14 (100%)	0 (0%)	0 (0%)	0 (0%)
04-Feb-20	53	53 (100%)	0 (0%)	0 (0%)	0 (0%)
09-Sep-20	4417	4400 (99.62%)	14 (0.32%)	2 (0.05%)	1 (0.02%)
12-Oct-20	5139	5126 (99.96%)	11 (0.21%)	0 (0%)	2 (0.04%)
05-Nov-20	3494	3478 (99.54%)	16 (0.46%)	0 (0%)	0 (0%)
04-Dec-20	3407	3397 (99.71%)	10 (0.29%)	0 (0%)	0 (0%)
04-Jan-21	3540	3530 (99.72%)	10 (0.28%)	0 (0%)	0 (0%)
04-Feb-21	3962	3944 (99.55%)	18 (0.45%)	0 (0%)	0 (0%)
01-Mar-21	4024	3999 (99.38%)	25 (0.62%)	0 (0%)	0 (0%)
07-Apr-21	4025	3991 (99.16%)	34 (0.84%)	0 (0%)	0 (0%)
06-May-21	3923	3870 (98.65%)	53 (1.35%)	0 (0%)	0 (0%)
01-Jun-21	3883	3820 (98.38%)	63 (1.62%)	0 (0%)	0 (0%)
06-Jul-21	3883	3823 (98.45%)	59 (1.52%)	0 (0%)	1 (0.03%)
02-Aug-21	3782	3723 (98.44%)	58 (1.53%)	0 (0%)	1 (0.03%)
01-Sep-21	3534	3476 (98.36%)	56 (1.58%)	2 (0.06%)	0 (0%)
01-Oct-21	3559	2882 (98.76%)	44 (1.24%)	0 (0%)	0 (0%)
01-Nov-21	3572	3527 (98.74%)	44 (1.23%)	1 (0.03%)	0 (0%)
06-Dec-21	3386	3527 (99.03%)	33 (0.97%)	0 (0%)	0 (0%)
03-Jan-22	3472	3442 (99.14%)	30 (0.86%)	0 (0%)	0 (0%)
01-Feb-22	3201	3170 (99.03%)	30 (0.94%)	0 (0%)	1 (0.03%)
01-Mar-22	3247	3222 (99.23%)	24 (0.74%)	0 (0%)	1 (0.03%)
01-Apr-22	2967	2947 (99.33%)	20 (0.67%)	0 (0%)	0 (0%)
01-May-22	2822	2802 (99.29%)	19 (0.67%)	1 (0.04%)	0 (0%)
01-Jun-22	2905	2883 (99.24%)	21 (0.72%)	1 (0.04%)	0 (0%)
01-Jul-22	2817	2793 (99.15%)	24 (0.85%)	0 (0.00%)	0 (0%)
01-Aug-22	2817	2793 (99.15%)	24 (0.85%)	0 (0.00%)	0 (0%)
01-Sep-22	2818	2794 (99.18%)	22 (0.78%)	1 (0.04%)	0 (0%)

As in any other primer used in qPCR techniques, the single mismatch caused by the Delta point mutation is not expected to prevent primer binding or functionality toward the SARS-CoV-2 genome. In theory, a single CoPrimer with 2 mismatches is expected to retain sensitivity with marked Ct delay and with 3+ mismatches, a CoPrimer is expected to have serious impairment (U.S. FDA, 2021) (FDA, 2021).

In addition, the qualitative detection of SARS-CoV-2 by more than 1 redundant marker minimizes the expected impact of a point mutation in any single CoPrimer or assay. As demonstrated by **Table 2**, only COVID-K-001 has a single marker, all the other products, COVID-K-002, ABC-K-001, and COVDS-K-003 / COVDS-K-004 target at least 2 markers.

As discussed previously, although the Co-Diagnostics risk analysis did not indicate that loss of sensitivity was likely with any of the products listed in **Table 2**. The impact of the single point mutation in lineage Delta was analyzed for T<sub>m</sub> impact on COVID-K-001. The impact on T<sub>m</sub> was determined to be modest, and the affected portion of the CoPrimer retained a predicted annealing temperature above that used in the validated thermocycling protocols. Therefore, the predicted T<sub>m</sub> impact analysis corroborated the previous risk determination that sensitivity toward Delta was unlikely to be affected. As a final analysis, synthetic RNAs were obtained with the Wild Type sequence and with the single point mutation present in the Delta lineage. The wet testing analysis confirmed that both RNAs exhibited estimated and confirmed Limits of Detection (LoD) within the 3-fold limit set as the acceptance criterion. Also, (King, et al., 2022) reportedly used the Logix Smart COVID-K-001 assay to accurately quantify COVID strains throughout the period from January 2021 through November 2021 as Delta became the predominant circulating strain.

Based on these evaluations, Co-Diagnostics remains confident that all Logix Smart SARS-CoV-2 assays, namely COVID-K-001, COVID-K-002, ABC-K-001, and COVDS-K-003/COVDS-K-004, retain full sensitivity for lineages classified as Variants of Concern (VOC) or Variants of Interest (VOI) at this time, including Omicron and all of its subvariants (B.1.1.529, B.1.1.529.XE, BA.1, BA.1.1, BA.2, BA.3, BA.4, BA.5, B.A.2 + L452X, B.A. 2.75), Delta, as well as other emerging variant sequences reflected in monthly analyses of the Nextstrain database thus far.

## References

- CDC. (2022, Feb 2). *COVID-19: SARS-CoV-2 Variant Classifications and definitions > Variant of Concern*. Retrieved Feb 2, 2022 from Centers for Disease Control and Prevention: <https://www.cdc.gov/coronavirus/2019-ncov/variants/variant-info.html>
- Faria, N. R., Claro, I., Candido, D., Franco, L., Andrade, P., Coletti, T., . . . Ge, o. b. (2021, Jan 12). *Genomic characterisation of an emergent SARS-CoV-2 lineage in Manaus: preliminary findings*, Sampled in 26-Jan-2021. Retrieved Jan 26, 2021 from Virological: <https://virological.org/t/genomic-characterisation-of-an-emergent-sars-cov-2-lineage-in-manaus-preliminary-findings/586>
- FDA. (2021, Jun 3). *SARS-CoV-2 Viral Mutations: Impact on COVID-19 Tests*. Retrieved Jun 24, 2021 from U.S. Food & Drug Administration: <https://www.fda.gov/medical-devices/coronavirus-covid-19-and-medical-devices/sars-cov-2-viral-mutations-impact-covid-19-tests>
- GISAID. (2020, Dec). *GISAID in the News: UK reports new variant, termed VUI 202012/01*. Retrieved Dec 21, 2020 from GISAID: <https://www.gisaid.org/references/gisaid-in-the-news/uk-reports-new-variant-termed-vui-20201201/>
- King, K. L., Wilson, S., Napolitano, J. M., Sell, K. J., Rennert, L., Parkinson, C. L., & Dean, D. (2022, May 10). SARS-CoV-2 variants of concern Alpha and Delta show increased viral load in saliva. *PLoS One*. doi:10.1371/journal.pone.0267750

- MHRA. (2022, Jan 28). Retrieved Feb 2, 2022 from <https://www.gov.uk/government/publications/covid-19-variants-genomically-confirmed-case-numbers/variants-distribution-of-case-data-28-january-2022>
- PANGO lineages. (2021, May 19). *Grinch | Global Report Investigating Novel Coronavirus Haplotypes: B.1.167.2*. Retrieved Jun 24, 2021 from PANGO lineages: [https://cov-lineages.org/global\\_report\\_B.1.617.2.html](https://cov-lineages.org/global_report_B.1.617.2.html)
- Public Health England. (2020, Dec 14). *PHE investigating a novel variant of COVID-19: A new variant of the virus that causes COVID-19 (SARS-CoV-2) has been identified across the South East of England*. Retrieved Dec 21, 2020 from GOC.UK: <https://www.gov.uk/government/news/phe-investigating-a-novel-variant-of-covid-19>
- Rambaut, A., Loman, N., Pybus, O., Barclay, W., Barret, J., Carabeli, A., . . . Volz, E. o.-1.-U. (2020, Dec 18). *Preliminary genomic characterisation of an emergent SARS-CoV-2 lineage in the UK defined by a novel set of spike mutations*. Retrieved Dec 21, 2020 from Virological: <https://virological.org/t/preliminary-genomic-characterisation-of-an-emergent-sars-cov-2-lineage-in-the-uk-defined-by-a-novel-set-of-spike-mutations/563/1>
- U.S. FDA. (2021, Jan 8). *Genetic Variants of SARS-CoV-2 May Lead to False Negative Results with Molecular Tests for Detection of SARS-CoV-2 - Letter to Clinical Laboratory Staff and Health Care Providers*. Retrieved Feb 2, 2021 from U.S. Food & Drug Administration: <https://www.fda.gov/medical-devices/letters-health-care-providers/genetic-variants-sars-cov-2-may-lead-false-negative-results-molecular-tests-detection-sars-cov-2>
- Virological. (2020, Dec 18). *Preliminary genomic characterisation of an emergent SARS-CoV-2 lineage in the UK defined by a novel set of spike mutations*. Retrieved Dec 21, 2020 from Virological: <https://virological.org/t/preliminary-genomic-characterisation-of-an-emergent-sars-cov-2-lineage-in-the-uk-defined-by-a-novel-set-of-spike-mutations/563>
- WHO. (2022, Jan 7). *Enhancing Readiness for Omicron (B.1.1.529): Technical Brief and Priority Actions for Member States Update #5*.
- WHO. (2022, Feb 2). *Tracking SARS-CoV-2 Variants*. Retrieved Feb 2, 2022 from World Health Organization: <https://www.who.int/en/activities/tracking-SARS-CoV-2-variants/>

In case of additional questions, please contact [support@codiagnostics.com](mailto:support@codiagnostics.com).



Kenneth K.C. Bramwell, Ph.D. Deputy Chief Science Officer

[k.bramwell@co-dx.com](mailto:k.bramwell@co-dx.com)

Co-Diagnostics, Inc.

2401 S. Foothill Dr, Ste. D, SLC, Utah 84109

[co-dx.com](http://co-dx.com)